

Release Notes for Version 1.0.0-Beta.0 of EmberEquine for Windows 10

February 2018

These notes provide important information for Juniper Systems release of the EmberEquine operating on any device running Microsoft Windows 10. These notes may also include important hardware information for the device.

Documentation and software updates are available from:

www.junipersys.com/emberequine/support

Known Issues

These issues are open in this version of the software.

- EmberEquine can freeze or lock when setting a custom Wi-Fi access point name to the same name as an existing Wi-Fi access point (297)
 - Do not use the skip button when completing multi-entry fields (such as L Ovary, R Ovary, L Cyst, & R Cyst) in the exam forms because it may cause EmberEquine to crash when viewing the exam history screen or when running reports (291)
 - Reports do not show results on a secondary contact (300)
 - Events that are not linked to an Animal are not found in an event report (171)
 - Having more than one user logged into Windows and in EmberEquine on the same PC is not supported (286, & 282)
-

PN 27689-00. © Copyright 2018. Juniper Systems, Inc. All rights reserved.

Windows is a registered trademark of Microsoft, Inc. in the United States and/or other countries.

Other trademarks in this document are the property of their respective owners.